

Reagent Provided

For Research Use Only. Not for use in diagnostic procedures.

Probe: 110 µl per vial. Ready to use.

The probe is premixed in hybridization buffer (formamide, dextran sulphate, SDS and SSC).

Warning and Precaution: Handle all reagents and slides containing fluorophores in reduced light to prevent photobleaching. The probe mixture contains formamide, avoid exposure. Handle all reagents with care; wear appropriate personal protective equipment.

Reagents & Material Required but Not Provided

Coplin jar	1N Hydrochloric acid (HCl)
Waterbath	Pepsin
Microcentrifuge	10X Phosphate Buffer Saline (PBS)
Centrifuge tube (15 ml)	20X Saline Sodium Citrate (SSC)
Slide tray	10% Formalin
Rubber cement	100% Ethanol
Cover glass	Tween 20
Forceps	DAPI/ Antifade
Temp. controlled hot plate	Autoclaved distilled water
Micropipettes (1-200 µl)	Acetic acid
MgCl ₂	Methanol

Reagent Preparation

2X SSC. Mix 100 ml of 20X SSC and 900 ml of distilled water. Store at room temperature (RT).

2X SSC/0.1% Tween 20. Add 100 ml of 20X SSC and 1 ml of Tween 20 to 900 ml distilled water. Mix well by swirling. Store at RT.

0.5X SSC/0.1% Tween 20. Add 25 ml of 20X SSC and 1 ml of Tween 20 to 975 ml distilled water. Mix well by swirling. Store at RT.

Ethanol Series (70%, 85%, and 100%). Prepare v/v dilutions of 100% ethanol with distilled water. Store at RT.

0.01N HCl. Add 0.5 ml of 1N HCl to 49.5 ml distilled water.

0.4% Pepsin Stock Solution (4 mg/ml). Dissolve 250 mg of pepsin in 62.5 ml of 0.01 N HCl. Store 2 ml aliquots at -20°C.

1X PBS. Mix 100 ml of 10X PBS to 900 ml distilled water. Adjust pH to 7.0. Store at RT.

Formalin Fixative. Add 12.5 ml of 10% formalin to 37 ml of 1X PBS and 0.5 ml of 100X MgCl₂. Store at RT.

3:1 Fixative. Add 15 ml of methanol to 5 ml acetic acid. Make fresh each time, discard unused portion.

Filter Requirements

Fluorophore	Excitation _{max}	Emission _{max}
DAPI	360 nm	460 nm
Aqua	431 nm	480 nm
Green	496 nm	551 nm
Gold	525 nm	551 nm
Red	580 nm	603 nm

Storage conditions: Hybridized slides and Probe must be stored at -20°C protected from direct light.

Protocol

Ready to use product.

Do not reconstitute or dilute with hybridization buffer.

Slide Preparation

1. Transfer the cervical cells stored in liquid-based medium to a 15 ml centrifuge tube and centrifuge at 1,200 rpm for 10 min at RT.
2. Discard the supernatant and resuspend the pellet in freshly made 3:1 fixative.
3. Repeat the wash process (centrifuge-resuspend in fixative), about 2-3 times to clear the cell pellet of all debris. Resuspend in 2 ml of fixative.

Procedure Note: Some specimens may require more or less fixative depending on the pellet size.

4. Prepare the slide according to the standard laboratory protocol. If not hybridizing within a few days, store the slide at 4°C, or -20°C for long-term storage (over a month).

Slide Pretreatment with Pepsin

1. Dilute the pepsin stock to a final concentration of 0.5 mg/ml using 0.01 N HCl. Incubate the slide for at least 10 min at 37°C.

Procedure Note: Some specimens may require longer digestion time in pepsin.

2. Wash the slide 5 min in 1X PBS at RT.
3. Incubate the slide for 5 min in 1% formaldehyde fixative at RT.
4. Wash the slide for 5 min in 1X PBS at RT.
5. Dehydrate the slide in 70%, 85%, and 100% ethanol at RT for 1 min each. Air dry.

Procedure Note: Check the morphology of the sample under phase contrast microscope before hybridization.

Probe Denaturation / Hybridization

1. Vortex the Probe briefly and spin the tube in a microcentrifuge.
 2. Apply 10 µl of the Probe, cover with a cover glass (22 x 22 mm).
- Procedure Note: Care should be taken to avoid air bubbles. Smaller or larger cover glass may be used with proportional change in Probe volume.*
3. Seal the edges of the cover glass thoroughly with rubber cement.
 4. Co-denature the slide and the Probe mixture for 3 min at 80°C on a temperature controlled hot plate protected from direct light.
 5. Incubate for 48 hours in a humidified environment at 37°C and protected from light.

Post Hybridization Washing

Procedure Note: Do not allow the slide to dry before the washes are complete.

1. Remove the rubber cement with forceps.
2. Remove the cover glass by soaking the slide in 2X SSC at RT.
3. Wash the slide 2 x 5 min in 2X SSC/0.1% Tween 20 at 45°C with agitation.
4. Briefly rinse the slide in distilled water and air dry the slide out of direct light.
5. Apply 20 µl of DAPI/ Antifade solution to the hybridized area and cover with a cover glass (25 x 25mm).

Note: Dilute the DAPI to 1:10 with Antifade to prevent the counterstain from overwhelming the intensity of the blue signal.

Contact Us:

Cancer Genetics Italia S.r.l.
Viale Luigi Majno, 17
20122 Milano – Italia

e-mail: support@cancerogeneticsitalia.com
www.cancerogeneticsitalia.com