

**ExpressArt Mag FFPE Clear RNAREady kit**

**FFPE RNA isolation kit**

**Catalogue No. 9011-A100**

**for 100 RNA isolations from archival paraffin samples**

**with Magnetic Beads for RNA isolation**

**No separate paraffin removal step !**

**This protocol provides the required laboratory procedures**

**For additional information about ExpressArt technology and products  
see**

**<http://www.amp-tec.com/products>  
» more information**

## Kit Contents

Box I	Volume per RNA isolation	Kit contents for 100 isolations
Lysis Buffer (LB)	870 µl	87 ml
FFPE Clear a formulation of <i>non-toxic, non-odorous</i> organic solvents	500 µl	55 ml
Demodification Solution (DeS)	130 µl	15 ml
RNase-free water	50 µl	5 ml

Box II	Volume per RNA isolation	Kit contents for 100 isolations
NucleoGuard	5 µl	0.6 ml
Proteinase K, PCR Grade*	30 µl	powder
RNase-Free DNase**	45 µl (working solution)	powder
DiB, Dilution buffer for DNase	39.13 µl	2x 2 ml
Magnetic Beads Stock	100 µl	10 ml

\* Dissolve proteinase K (powder) in 4 ml water (25 mg/ml); store at -20°C.

\*\*Prepare “DNase working solution”

**First step:**

Dissolve DNase (powder) in 600 µl water; store at -20°C.

**Second step:**

**DNase working solution** (prepare according to your needs)

The example is for 51 isolations:

Mix 0.3 ml DNase + 2.0 ml Dilution Buffer (DiB)

Unused Aliquots are stored at -20°C

**Please note, these reagents are not included in the kit**

Isopropanol

### Required instruments

Thermomixer (Eppendorf #5350 000.013 or “Thriller” from PEQLAB #91-7010) or equivalent

Magnetic separator

### Tissue samples

**Please note:**

1) to avoid negative surface effects (contact to air), do not use the first 2 slides for preparation of nucleic acids.

2) no need for immediate processing of slides,


freezing and short-term storage (up to 3 days) at -80°C is even advantageous:

there is no negative effect on RNA quality; frequently, RNA yields will increase significantly

Up to 5 FFPE slides (10 µm thick with appr. 1-2 cm<sup>2</sup> of tissue) per isolation can be used in the lysis procedure

### Work Flow

#### RNA isolation from FFPE samples using the ExpressArt Mag FFPE Clear RNAREady kit



**Total Time: ~2 h**

### **1. Sample preparation**

Insert paraffin sections in reaction tube

Immediately (**no prior paraffin removal required**) add 500 µl **FFPE Clear**

**(Note:** For immediate processing, it can be more convenient to prepare tubes with FFPE Clear, and transfer paraffin sections directly into the solvent, avoiding "electrostatic problems")

**Heat for 1 min @ 80°C.**

Add 500 µl **Lysis Buffer (LB)** + 30 µl Proteinase K (25 mg/ml) + 5 µl NucleoGuard

**2. Lysis:** Incubation in Thermomixer (1000 rpm) for 1 h @ 55°C

**3. Demodification:** Add 130 µl **Demodification Solution (DeS)** (to 535 µl lysate)

Vortex, then incubation in Thermomixer (1000 rpm) for 20 min @ 80°C

Chill on ice (3 min)

### **4. RNA Purification**

**Before you start, prepare "bead suspension"**

**For 1 sample:** Transfer 100 µl of Magnetic Beads Stock into 1.5 ml tube, place on magnet, discard liquid and resuspend pellet in 900 µl 70% isopropanol.

**For multiple samples,** use two suspension steps

Example for **10 samples:** **Step 1:** Transfer 1 ml of Magnetic Beads Stock into 1.5 ml tube, place on magnet, discard liquid and resuspend pellet in 1 ml of 70% isopropanol by pipetting.

**Step 2:** Transfer into 15 ml tube, add 8 ml 70% isopropanol and mix thoroughly.

#### **4.1. Recovery of the aqueous lysate**

Centrifuge lysate (from step 3): 3 min @ 13'000 rpm

This results in three zones: pellet at the bottom (debris), **aqueous phase** with the lysate in the middle and organic supernatant at the top.

#### **4.2. RNA binding to beads**

Transfer 450 µl of the **aqueous phase** (avoid transfer of the organic supernatant) in 1.5 ml reaction tube

Mix 450 µl of the aqueous lysate with 900 µl bead suspension by pipetting (= 1350 µl)

Incubate on Thermomixer (750 rpm) for 5 min @ room temperature

Transfer on magnet, discard supernatant

#### **4.3. First washing steps**

Wash 2x with 800 µl 70% Isopropanol: Mix by pipetting, place on magnet, remove liquid

**For optional DNA-Digestion** proceed with steps 4.4 to 4.6, or go directly to step 4.7.

#### **4.4. DNA-Digestion**

Incubate on Thermomixer (900 rpm) for 2 min @ room temperature to dry the beads.

Add 45 µl of the "DNase Working solution" to beads, mix by pipetting.

Incubate on Thermomixer (750 rpm) for 15 min @ room temperature

#### **4.5. DNase inactivation and RNA recovery**

Add 355 µl Lysis Buffer + 800 µl 70% Isopropanol and mix by pipetting

Incubate on Thermomixer (750 rpm) for 5 min @ room temperature

Transfer on magnet, discard supernatant

#### **4.6. Second washing steps**

Wash 2x with 800 µl 70% Isopropanol: Mix by pipetting, place on magnet, remove liquid

#### **4.7. Elution**

Incubate on Thermomixer (900 rpm) for 3 min at 65°C to dry the beads completely

Add 50 µl water, mix by pipetting and make sure to soak the beads completely

Incubate on Thermomixer (900 rpm) for 5 min @ 65°C

Transfer on magnet and transfer supernatant in fresh tube

**RNA is recovered in ~45 µl.**

**Expected yield 2 - 6 µg RNA per FFPE slide**

#### **AmpTec GmbH**

Koenigstrasse 4a | 22767 Hamburg | Germany

Tel: +49 (0)40 – 636 747 22 | Fax: +49 (0)40 – 636 747 19

Technical Support:

"Dr. Guido Krupp" <krupp@amp-tec.com>

"Dr. Peter Scheinert" <scheinert@amp-tec.com>

www.amp-tec.com